





TA74 4G T



The TA74 product range offers multiple innovations that help you manage your telecare service easily, efficiently and economically. The devices are equipped with an internal 4G module which enable communication via the mobile network.

Two "plug and play" programming methods are now available, one through the TeleAlarm® Cloud Services platform and the other through the Configuration Manager software.

System overview



- 1. Daily button
- 2. Action button
- 3. Sign in / sign out button
- 4. Emergency call button
- 5. Maintenance button

Functions

Multi-purpose buttons

In case of an emergency, the TA74 sends an alarm via the telephone network, mobile network or via LAN/Ethernet connection, either to an alarm receiving centre or a private phone connection.

The emergency call is triggered by pressing the emergency call button on the TA74, or on the button of the wireless transmitter that is delivered with the device. Up to 20 wireless transmitters can be paired to the TA74. The action button is a multi-purpose trigger: It can be pressed to cancel an unintentionally triggered alarm during the pre-

- Quick and easy programming via intuitive cloud programming
- Various communication protocols
- Vocal guidance
- First-class audio quality
- Best of modern communication options with Ethernet/VoIP, PSTN and GSM/GPRS/UMTS/LTE/ VoLTE
- Always up to date via remote software updates
- Activity monitoring
- Vocal reminders
- Easy care easy to clean and waterproof in accordance with IP32
- 100% high-quality ABS material in premium white

alarm time, but it can also be configured to be used as a service button.

When the TA74 has an error (message), pressing the action button repeats the status announcement.

Pressing the sign in / sign out button tells the TA74 whether the user is currently at home or not. Finally, the TA74 includes a built-in activity monitor. This ensures that an emergency call is sent automatically

- if the Daily button is not pressed,
- if no activity is detected during a preset period of time.

Multiple protocols

The TA74 is capable of multiple protocols and therefore compatible with most telecare monitoring solutions.

Multiple programming possibilities

The TA74 can be programmed via the TeleAlarm[®] Cloud Services platform (over LAN/Ethernet, analogue telephony or mobile network)) or by Configuration Manager software on a local computer.

This programming method is a "plug and play" process, sparing installers time and effort.

The TA74 can also be programmed by compatible alarm receiving centres.

Vocal reminders

Three reminder messages can be recorded.

Each notification can be automatically played on a configured day and at a configured time. It is repeated every 20 seconds until the Daily button is pressed or the preset time has elapsed. If the reminder message is not acknowledged, the device can send an alert to the alarm receiving centre.

Interactivity with wireless transmitters and detectors

A wide range of wireless transmitters and detectors can be used with the TA74, for example the Tx Transmitter, the wireless motion detector or the wireless smoke detector.





TA74 technical specifications

Discussion	F0 101 101
Dimensions (H × W × D)	50 × 191 × 191 mm
Weight	Approx. 0.6 kg
Permissible temperature range (environment / storage)	+5 °C to +40 °C –20 °C to +60 °C
Maximum relative humidity	93 %
Altitude	Up to 5000 m
Switched-mode power supply	AC adaptor 100 - 240 V AC, 50/60 Hz, 0.4 A Input, 6.1 V DC, 1.6 A Output
Current consumption in standby mode	Approx. 150 mA with power supply. Approx. 38 mA on battery
Unit battery	3.6 V NiMH battery, 2.2 Ah
Battery operating time	Min. 36 h with one 30-minute call. All values on date of purchase with a fully charged battery
Communication connections	Ethernet/VoIP, PSTN, GSM/ GPRS/UMTS/LTE/VoLTE
Receiver	Complies with category 1 as specified by EN 300220-1 V3.1.1
Frequency	869.2125 MHz
Environmental class	Class 1 as specified by EN 50134- 3
Degree of protection	IP32 (IP30 wall mounting)
Max. number of wireless transmitters and detectors	20
Programmable call destinations	10 call destinations freely assignable to any type of trigger group
Voice message	Recordable, max. 10 seconds
Volume settings	8 levels
Emergency call protocols	Multi-protocol (TS 50134-9, RBIP-S, SCAIP, RBIP, SIP, telephone, RB2000, RB2000E, TTnew+, CPC, BS8521)
Programming options	 Remote programming via TeleAlarm[®] Cloud Services Configuration Manager Remote programming from alarm receiving centre Service mode
Hands-free range	Min. 15 m (voice possibility when telephone is on-hook)
Self-test function/ test call	 Analogue protocols: Freely selectable from 1 to 28 days. IP protocols: Freely selectable from 1 to 999 hours
Available languages	Dutch, English, French, German, Spanish, and tone only

User guidance	Synthetic speech and LED indication
User assistant	"Synthetic speech disabled" mode
Connectivity	 LAN/Ethernet connection (10/100 Mbps) Phone line (PSTN) External telephone micro-SIM card socket
SIM card size	• Micro-SIM (15 mm × 12 mm)
Mobile network	2G GSM: GSM900, DCS1800 3G UMTS: B1 (1200), B8 (900) 4G LTE: B1, B3, B7, B8, B20, B28
Vocal reminders	Three recordable reminder messages (max. 15 seconds) (programming of day/time of playback)

Technical specifications of wireless transmitter

Dimensions (H × D)	11 × 41 mm
Weight	15 g
Permissible temperature range	Operation: -10 °C to +55 °C Storage: -20 °C to +80 °C
Electrical power supply	Replaceable lithium battery CR2430
Current consumption	Transmission power: max. 22 mA Quiescent mode: max. 3 µA
Battery lifetime	Approx. 5 years if one emergency call is triggered per day. Automatic radio link and battery monitoring every 21 hours
Frequency	869.2125 MHz
Stability (0 °C to +55 °C)	< 2.5 kHz
Modulation	FSK
Antenna	Integrated
Radiated intensity	Approx. 320 µW
Degree of protection	IP67 in accordance with IEC 60529
Environmental class	II
Maximum range	 Inside buildings: 30 – 50 m (depending on the building structure) In the open air: max. 300 m





Certifications and approvals

LVD 2014-35-EU, EMC 2014-30-EU, RED 2014/53/EU, RoHS 2011/65/EU. EN 50134-1, EN 50134-2, EN 50134-3, EN 50134-5 EN 50134-9 EN 301489-1 V2.1.1 EN 301489-3 V2.1.1 EN 62368-1 EN 62479 EN 50130-4 EN 55032 EN 301 511 V12.5.1 EN 301 908-2 V11.1.1 EN 301 489-7 V1.3.1 EN 301 489-24 V1.5.1 EN 300220-2 V3.1.1 Category 1 radio receiver Conducted sensitivity typ. –115 dBm Field sensitivity < 36 dBµV/m

Scope of delivery

Quantity	Component
1	Pre-installed rechargeable battery
1	Built-in 4G module: 2G dual band 900/1800 MHz, UMTS (3G) B1/B8 and LTE (4G) B1/B3/B7/B8/B20/B28
1	Wireless transmitter including CR2430 battery and necklace or bracelet
1	Power adapter
1	Optional telephone cable and/or network connection cable
1	Optional manual for technicians/installers
1	Ouick Guide for users
1	Optional manual for technicians/installers
1	Quick Guide for users

TeleAlarm SA Rue du Pont 23 2300 La Chaux-de-Fonds Switzerland www.telealarm.com

© Copyright TeleAlarm SA, 2023

TeleAlarm Europe GmbH Hertzstraße 2 04329 Leipzig Germany

Contact in UK: Tel.: +44 (0) 333 0124392 info-uk@telealarm.com