

Title	Status	Date	Version
L071RX FW description file	Final	17.10.2022	1.2

1 Forewords

This document applies to the NurseCall L071RX FW programming, dedicated to the LE80, AP80 and NC8 products.

Files are supplied in binary format.

2 L071RX firmware:

Software SAP#:	T.040.001.513
Software SAP Description:	Software LE80
Software information to be printed	None

2.1 Binary firmware version details:

File format:	binary
Version:	NC_05.A00.02.1476
Build date:	17.10.2022
Compiler version	IAR EWARM 8.50.6
File Name:	NC_05.A00.02.1476.bin
MD5 Checksum:	35DB3BD1DD22B9E12FBF28A9F2A3858F
Remark:	None

2.2 Change History:

SW Version	Build date	Description
05.A00.01.1410	13.07.2021	Ready for Ramp-up production
05.A00.01.1422	01.11.2021	Fixed character problem in Production initialization
05.A00.02.1476	17.10.2022	NC8: ACK by B81/B80A generates ACKN Patient + ACKN Badge instead of ACKN Patient + ALARM Badge. AP80: 1) Possibility to disable the forwarding of the "Supervisions" to the serial output, avoiding too much information. 2) Filter RF-FIFO Buffer to avoid sending several same serial datagrams when First level of RF retries.